

Support and Training



We offer full certificated training for all your operational team to ensure they are confident and competent to use our machines.



Our friendly and knowledgeable customer care and engineering teams are on hand to support you when you need it.



A full range of maintenance validation packages are available to be tailored to the needs of your organisation.

Let's raise the bar together -
aiming for higher industry standards

K+F machines are compliant with the MDD and MDR, along with the standard BS EN ISO 15883 Parts 1, 5 & 7. BES are also providing a decontamination health check and road map to compliance for your service.

For more information, a quote, or to book a free site survey, get in touch with us on 0117 966 67 61 or email info@beshealthcare.net

K+F Equipment Washer-Disinfectors

Compliant with BS EN ISO 15883 Parts 1, 5 & 7

Now with
integrated
dryers as
standard



KF900



B | E | S
HEALTHCARE

0117 966 6761
info@beshealthcare.net
www.beshealthcare.net



B | E | S
HEALTHCARE

0117 966 6761
info@beshealthcare.net
www.beshealthcare.net



Effortless Cleaning and Disinfection

K+F offers a range of washer-disinfectors designed specifically for community equipment, providing an automated, cost-effective solution for cleaning and disinfecting a variety of medical devices such as wheelchairs, hoists, bed frames, and commodes. Each machine includes an integrated drying facility for added convenience.

With extensive expertise in the design, production, and manufacturing of cleaning and disinfection systems, K+F supports hospitals, community equipment services, and decontamination departments with reliable, high-quality solutions tailored to their needs.



KF300



Compliant Machines

Complies with the Medical Device Regulation class IIa, meeting BS EN ISO 15883 Parts 1,5 and 7.



Cost Saving

Saves costs in relation to time spent, energy, staffing, equipment replacement in comparison to manual washing.



Time Saving

Equipment can be processed rapidly compared to the longer time it takes to undertake manual cleaning.



Peace of Mind

Effective, repeatable cleaning and disinfecting cycles ensure your equipment is reliably safe to reissue.

The K+F Range

KF200
WASHER
DISINFECTOR



Ideal for washing wheelchairs, hoists, and other items such as walking frames.

KF300
PASS THROUGH
WASHER



Designed to wash and disinfect foldable beds, and larger pieces of equipment.

KF900
PASS THROUGH
WASHER



Developed to wash and disinfect beds. Ideal for those with a high level of throughput.

	KF200	KF300	KF900
External Length (cm)	138	138	250
External Width (cm)	164	164	211
External Height (cm)	200	250	250
Chamber Length (cm)	135	135	245
Chamber Width (cm)	90	90	143
Chamber Height (cm)	140	190	190