



# chunc one

Manual Attendant Wheelchair For Children & Adults

## Key Features

- Excellent all-rounder
- Highly adaptable
- Indoor & outdoor wheelchair
- Long-term solution
- Available in five sizes: XS to XL
- Up to 110kg user weight
- Crash tested to ISO 7176-19

## What's Included

- Tilt in space chair
- Push handles
- Foot operated parking brakes
- Foot operated anti-tip levers

## Optional Extras

- Recline to 180-degrees
- Hand-operated drum brakes
- Postural support
- Medical & lifestyle accessories



**BES**  
HEALTHCARE  
IRELAND

+44 (0) 1179 666 761  
+353 (0) 216 019 950  
[info@beshealthcare.ie](mailto:info@beshealthcare.ie)  
[www.beshealthcare.ie](http://www.beshealthcare.ie)



## Uniquely adaptable, tailored to your needs

Chunc One is the signature chair from Chunc. Durable, functional and stylish, it is suitable for those with moderate to complex disabilities.

Every chair can be adapted to suit the child's individual needs. Our seating specialist will work with your therapist on a one-to-one basis to understand your specific requirements.

## Grows with you

The Chunc One's design means that using our Growth Kits, your Chunc One can grow as your child grows or their needs change. In other words, you can keep one chair for life.

## Tilt-in-space wheelchair

Chunc One has a 45-degree tilt-in-space feature which means the chair can be tilted whilst maintaining correct posture. Being able to tilt the chair provides pressure relief as well as allowing a greater feeling of comfort and relaxation, when required.

## Easy to recline

The recline option enables the Chunc range of wheelchairs to be reclined back to a horizontal position – ideal when you wish to use the chair as a changing station, or even in a medical emergency.

SPECIFICATION	XS	SMALL	MEDIUM	LARGE	XL
Max User Mass	110kg	110kg	110kg	110kg	110kg
Seat Cushion Depth Range	30-35cm	35-40cm	40-45cm	45-50cm	50-55cm
Seat Width	30cm	35cm	40cm	45cm	50cm
Back Rest Height	40cm	45cm	51cm	57cm	62cm
Total Backrest Height Growth	5cm	5cm	5cm	5cm	5cm
Backrest Width	25cm	30cm	31.5cm	35cm	48cm
Leg Rest Angle Range	80°-170°	80°-170°	80°-170°	80°-170°	80°-170°
Leg Rest Length Range	21-48cm	21-48cm	21-48cm	21-48cm	21-48cm
Back Rest Recline Angle Range	90-180°	90-180°	90-180°	90-180°	90-180°
Tilt in Space	+3° to -45°	+3° to -45°	+3° to -45°	+3° to -45°	+3° to -45°



**BES**  
HEALTHCARE  
IRELAND

+44 (0) 1179 666 761  
+353 (0) 216 019 950  
info@beshealthcare.ie  
www.beshealthcare.ie

