



Height adjustable indoor OR outdoor wheelchair

Key Features

- Easy to adjust height, from floor to highchair level
- Perfect for school, home or for outdoor use
- Grows with you

What's Included

- Tilt in space chair
- Push handles
- Outdoor: Drum brake with parking lever, anti-tip levers, tie down points
- Indoor: Castors with built-in parking brake function

Optional Features

- Recline wheelchair to 180-degrees
- Postural support, medical & lifestyle accessories

INDOOR



OUTDOOR



(Colours for indoor & outdoor versions)



BES
HEALTHCARE
IRELAND

+44 (0) 1179 666 761
+353 (0) 216 019 950
info@beshealthcare.ie
www.beshealthcare.ie



Outdoor Hilo: Join in with every activities

The Outdoor Hilo ensures you can take part in a wide range of activities at home, at school and outside. It is quick and simple to adjust the height, using the remote control. This makes it easy to get involved socially at all levels, from working at a desk in school, to playing games or enjoying crafts at the table to making a cake at the kitchen unit.

The Outdoor Hilo’s chunky rear wheels mean users can also be included in outdoor activities, from school trips and family excursions to picnics and bbqs with friends.

Indoor Hilo: Get involved at every level

The Indoor Hilo makes it easy and comfortable to participate in a wide range of activities at home and in school. Whatever is happening, you can join in!

Whether you pick the foot pump or powered wheelchair, it is easy and quick to adjust the height of the chair. And thanks to its castor wheels, the Indoor Hilo is also simple to move from room to room.



	XS	SMALL	MEDIUM	LARGE
Max User Mass	75kg	75kg	75kg	75kg
Seat Cushion Depth Range	30 – 35cm	35 – 40cm	40 – 45cm	45-50cm
Seat Width	30cm	35cm	40cm	45cm
Back Rest Height	40cm	45cm	51cm	57cm
Total Back Rest Height Growth	5cm	5cm	5cm	5cm
Back Rest Width	25cm	30cm	31.5cm	35cm
Leg Rest Angle Range	80°-170°	80°-170°	80°-170°	80°-170°
Leg Rest Length Range	21-48cm	21-48cm	21-48cm	21-48cm
Back Rest Recline Angle Range	90-180°	90-180°	90-180°	90-180°
Tilt in Space	+3° to -35°	+3° to -35°	+3° to -35°	+3° to -35°



BES
HEALTHCARE
IRELAND

+44 (0) 1179 666 761
+353 (0) 216 019 950
info@beshealthcare.ie
www.beshealthcare.ie

