



# SPICA

## Rehabilitation Spica Manual Wheelchair

### Key Features

- Aids post-operative recovery
- Support for fixed cast positions
- Available in four sizes: XS to L
- Highly adaptable
- Indoor & outdoor wheelchair
- Up to 110kg user weight

### What's Included

- Tilt-in-space chair
- Removable, adjustable leg splints
- Foot operated parking brakes
- Foot operated Anti-tip levers
- Tie down points for safe travel

### Optional Extras

- Range of seat & backrests
- Recline to 180 degrees
- Postural Support options



(Only available in black)



**BES**  
HEALTHCARE  
IRELAND

+44 (0) 1179 666 761  
+353 (0) 216 019 950  
[info@beshealthcare.ie](mailto:info@beshealthcare.ie)  
[www.beshealthcare.ie](http://www.beshealthcare.ie)





## Aids recovery from spica surgery

Recommended by major NHS and private children's hospitals in the UK, the Spica is Chunc's post-operative wheelchair. Built primarily to accommodate hip spica casts, it facilitates early release from hospital after surgery to the familiar surroundings of home, whilst also freeing up hospital beds and saving the NHS money.

The Chunc Spica wheelchair is very flexible, with options to accommodate a wide range of fixed cast positions, from back recline and leg angle to knee bend. Angle adjustable footrests can also be fitted to assist physiotherapy. The Chunc Spica also suits those with other requirements, such as radial leg extension braces and spinal surgery.



## Tilt-in-space chair

The Chunc Spica comes with a 45-degree tilt in space feature for added comfort as well as to improve posture and hip position.

## Rental model available

Subject to quotation a rental model of the Spica is available. Contact our sales team for more information.

SPECIFICATION	XS	SMALL	MEDIUM	LARGE
Max User Mass	110kg	110kg	110kg	110kg
Seat Cushion Depth Range	30-35cm	35-40cm	40-45cm	45-50cm
Seat Width	30cm	35cm	40cm	45cm
Back Rest Height	40cm	45cm	51cm	57cm
Total Backrest Height Growth	5cm	5cm	5cm	5cm
Backrest Width	25cm	30cm	31.5cm	35cm
Leg Rest Angle Range	80°-170°	80°-170°	80°-170°	80°-170°
Leg Rest Length Range	21-48cm	21-48cm	21-48cm	21-48cm
Back Rest Recline Angle Range	90-180°	90-180°	90-180°	90-180°
Tilt in Space	+3° to -45°	+3° to -45°	+3° to -45°	+3° to -45°



**BES**  
HEALTHCARE  
IRELAND

+44 (0) 1179 666 761  
+353 (0) 216 019 950  
info@beshealthcare.ie  
www.beshealthcare.ie

