

EPIC SEATING ORDER FORM

Posterior Pelvic Support Sizes.

Width Options

Small: Fits 10"-16" frame width
Medium: Fits 16"-21" frame width
Large: Fits 21"-26" frame width

Height

The posterior pelvic support is 6 inches tall. The pad is adjustable along the cane, allowing the pad to rest between the iliac crest and the PSIS.

RESISTANCE LEVELS > LOW

Low to Normal Tone

Ideal for users with limited active movement who may also experience posterior pelvic tilt. Used to allow relaxation of postural muscles while providing an articulating posterior thoracic support or to allow for manual adjustment and accommodation of limited active movement or kyphosis without the use of tools.

< HIGH RESISTANCE LEVELS

Low to High Tone

Ideal for users with normal to highly active movement, such as dystonic users who may exhibit an extension pattern of spasticity. Provides an articulating posterior thoracic support which can absorb energy exerted by the increased tone.

REA AZALEA

LOW

HIGH

EPiC Seating ASM Small 50N. [ESAT021601050.](#)

EPiC Seating ASM Medium 50N. [ESAT021602050.](#)

EPiC Seating ASM Large 50N. [ESAT021603050.](#)

EPiC Seating ASM Small 100N. [ESAT021601100.](#)

EPiC Seating ASM Medium 100N. [ESAT021602100.](#)

EPiC Seating ASM Large 100N. [ESAT021603100.](#)

OTTOBOCK DISCOVERY

LOW

HIGH

EPiC Seating ASM Small 50N. [ESAT021601050.](#)

EPiC Seating ASM Medium 50N. [ESAT021602050.](#)

EPiC Seating ASM Large 50N. [ESAT021603050.](#)

EPiC Seating ASM Small 100N. [ESAT021601100.](#)

EPiC Seating ASM Medium 100N. [ESAT021602100.](#)

EPiC Seating ASM Large 100N. [ESAT021603100.](#)

LEGGERO ENZO

LOW

HIGH

EPiC Seating ASM Small 50N. [ESAT040401050.](#)

EPiC Seating REM ASM Small 50N. [ESAT040404050.](#)

EPiC Seating ASM Medium 50N. [ESAT040402050.](#)

EPiC Seating REM ASM Med 50N. [ESAT040405050.](#)

EPiC Seating ASM Large 50N. [ESAT040403050.](#)

EPiC Seating REM ASM Large 50N. [ESAT040406050.](#)

EPiC Seating ASM Small 100N. [ESAT040401100.](#)

EPiC Seating REM ASM Small 100N. [ESAT040404100.](#)

EPiC Seating ASM Medium 100N. [ESAT040402100.](#)

EPiC Seating REM ASM Med 100N. [ESAT040405100.](#)

EPiC Seating ASM Large 100N. [ESAT040403100.](#)

EPiC Seating REM ASM Large 100N. [ESAT040406100.](#)

QUICKIE IRIS

LOW

HIGH

EPiC Seating ASM Small 50N. [ESAT090901050.](#)

EPiC Seating ASM Medium 50N. [ESAT090902050.](#)

EPiC Seating ASM Large 50N. [ESAT090903050.](#)

EPiC Seating ASM Small 100N. [ESAT090901100.](#)

EPiC Seating ASM Medium 100N. [ESAT090902100.](#)

EPiC Seating ASM Large 100N. [ESAT090903100.](#)

KI FOCUS CR

LOW

HIGH

EPiC Seating ASM Small 50N. [ESAT030501050.](#)

EPiC Seating ASM Medium 50N. [ESAT030502050.](#)

EPiC Seating ASM Large 50N. [ESAT030503050.](#)

EPiC Seating ASM Small 100N. [ESAT030501100.](#)

EPiC Seating ASM Medium 100N. [ESAT030502100.](#)

EPiC Seating ASM Large 100N. [ESAT030503100.](#)

INVACARE SOLARA 3G

LOW

HIGH

EPiC Seating ASM Small 50N. [ESAT021601050.](#)

EPiC Seating ASM Medium 50N. [ESAT021602050.](#)

EPiC Seating ASM Large 50N. [ESAT021603050.](#)

EPiC Seating ASM Small 100N. [ESAT021601100.](#)

EPiC Seating ASM Medium 100N. [ESAT021602100.](#)

EPiC Seating ASM Large 100N. [ESAT021603100.](#)

FORCE

EPIC SEATING REPLACEMENT PARTS

| | | | |
|---|--|--|--|
| Replacement Posterior Pelvic Block. | 10 wide, 6in Tall. <i>Part: ESR0008010.</i> | 11 wide, 6in Tall. <i>Part: ESR0008011.</i> | 12 wide, 6in Tall. <i>Part: ESR0008012.</i> |
| | 13 wide, 6in Tall. <i>Part: ESR0008013.</i> | 14 wide, 6in Tall. <i>Part: ESR0008014.</i> | 15 wide, 6in Tall. <i>Part: ESR0008015.</i> |
| | 16 wide, 6in Tall. <i>Part: ESR0008016.</i> | 17 wide, 6in Tall. <i>Part: ESR0008017.</i> | 18 wide, 6in Tall. <i>Part: ESR0008018.</i> |
| | 19 wide, 6in Tall. <i>Part: ESR0008019.</i> | 20 wide, 6in Tall. <i>Part: ESR0008020.</i> | 21 wide, 6in Tall. <i>Part: ESR0008021.</i> |
| | 22 wide, 6in Tall. <i>Part: ESR0008022.</i> | | |
| | Posterior Pelvic Support. | Lateral Pelvic Wedge Kit. <i>Part: ESR0006000.</i> | |
| Replacement Struts. | Single Gas Strut, 50N. <i>Part: ESR0003050.</i> | | Qty |
| | Single Gas Strut, 100N. <i>Part: ESR0003100.</i> | | Qty |
| | Cable Assembly. <i>Part: ESR0004000.</i> | | Qty |
| Pivot Parts. | Replacement Release Pin. <i>Part: ESR0001000.</i> | | Qty |
| | Replacement Bearings. <i>Part: ESR0002000.</i> | | Qty |
| Attachment Tube (pair). | Invacare Solara 3G, Attachment Tube. <i>Part: ESR020000.</i> | | |
| | Ki Focus CR, Attachment Tube. <i>Part: ESR0300000.</i> | | |
| | Leggero Enzo, Attachment Tube. <i>Part: ESR0400000.</i> | | |
| | Quickie Iris, Attachment Tube. <i>Part: ESR0900000.</i> | | |
| Foam Grip Handles. | Replacement Foam Grip Handles. <i>Part: ESR0005000.</i> | | Qty |

NOTES & SPECIAL REQUESTS